

## **Discussion on *Staphylococcus schleiferi***

### **Taxonomy**

The *Staphylococcus* genus is comprised of 49 species, including 25 subspecies, classified in the family Staphylococcaceae, order Bacillales, and phylum Firmicutes.<sup>1</sup> Coagulase detection has long been used to separate the more virulent *Staphylococcus aureus* group from species that are coagulase negative. Moreover, some zoonotic species have also been identified as coagulase-positive or variable such as *Staphylococcus schleiferi*. *Staphylococcus schleiferi* can be further divided into two subspecies, the coagulase-negative *S. schleiferi* subsp. *schleiferi* and coagulase-positive *S. schleiferi* subsp. *coagulans*.<sup>2</sup>

### **Identification**

Staphylococci are gram-positive cocci measuring 0.5 to 1.5 microns in diameter that appear as single cells, pairs, irregular clusters and sometimes short chains. The majority are facultatively anaerobic, non-motile and non-spore-forming as well as catalase-positive and able to grow in 10% NaCl.<sup>2</sup> Most *Staphylococcus* spp. will grow within 18 to 24 hours of incubation on routine media such as sheep blood agar under aerobic conditions. Colonies of *S. schleiferi* are medium to large in size, round, non-pigmented and are beta-hemolytic on blood agar. Several biochemical reactions can be used to differentiate *S. schleiferi* from other *Staphylococcus* species that cause human infections. *Staphylococcus schleiferi* is characteristically pyrrolidonyl arylamidase (PYR) positive, alkaline phosphatase positive, ornithine decarboxylase negative, and susceptible to novobiocin.<sup>1</sup> The two subspecies can be differentiated by the slide coagulase (clumping factor), tube coagulase (free coagulase) and urease reactions. *S. schleiferi* subsp. *schleiferi* is clumping factor positive, tube coagulase negative and urease negative whereas, *S. schleiferi* subsp. *coagulans* is clumping factor negative, tube coagulase positive and urease positive.<sup>2</sup> An international, multicenter study found that 25-75% of *S. schleiferi* can give false-positive latex agglutination tests used for the identification of *S. aureus*.<sup>3</sup> Due to  $\beta$ -hemolysis on blood agar and positivity by coagulase tests (eg, slide, tube and latex agglutination methods), misidentification of *S. schleiferi* as *S. aureus* can occur. A PYR test could be performed and if positive should raise suspicion that an isolate is not *S. aureus*. *S. schleiferi* subsp. *coagulans* or *Staphylococcus intermedius* group organisms fit the description of a coagulase-positive staphylococci that are PYR positive. These species can be differentiated based on sugar fermentation patterns.<sup>1</sup>

Scarce data are available on the performance of commercially available biochemical identification systems to identify *S. schleiferi*. Additional testing, such as sugar fermentation, for confirmation of identification is often required with commercial biochemical systems.<sup>2</sup> Matrix-assisted laser desorption/ionization time-of-flight mass spectrometry (MALDI-TOF MS) systems have been demonstrated to reliably identify *S. schleiferi* to the species level but is unable to differentiate the subspecies.<sup>4,5</sup> Increasing reports of *S. schleiferi* correlate with the widespread implementation of MALDI-TOF MS that has allowed clinical microbiology laboratories to better identify staphylococci species, including *S. schleiferi*.

### **Clinical Significance**

*Staphylococcus schleiferi* is an increasingly reported zoonotic pathogen that can cause opportunistic infections in humans. It colonizes the skin and mucosal surfaces of small animals, including dogs, cats, parrots and pigs.<sup>1,2,6</sup> In dogs, *S. schleiferi* causes canine otitis media, otitis externa and pyoderma.<sup>2</sup> As in this case, *S. schleiferi* infections in humans are often associated with wound infections from dog bites. Of the subspecies, *S. schleiferi* subsp. *schleiferi* is more commonly associated with human infections, including skin and soft tissue infections, device infections, osteomyelitis, and bacteremia. Few cases of human infection caused by *S. schleiferi* subsp. *coagulans* have been documented, although wound infections, endocarditis, and device infections have been reported.<sup>2,6</sup> Risk factors associated with infection include immunosuppression, malignancies, and recent surgical procedures suggesting nosocomial acquisition. However, the association of *S. schleiferi* with small animals and human infections supports the possibility that close contacts with these animals (pets or livestock) can result in transmission from animals to humans.<sup>1,2</sup>

## Antimicrobial Resistance and Therapy Considerations

From 39-73% of veterinary *S. schleiferi* isolates have been reported as oxacillin resistant.<sup>6</sup> Similar to other staphylococci, oxacillin resistance infers resistance to all beta-lactam antibiotics except for new anti-methicillin-resistant *S. aureus* cephalosporins (eg, ceftaroline). Beta-lactam resistance, is mediated by *mecA* (PBP2a) and has been associated with the staphylococcal cassette chromosome *mec* types I and IV in *S. schleiferi*.<sup>2,6</sup> In addition to oxacillin resistance, reduced susceptibility to erythromycin, clindamycin and fluoroquinolones has been reported.<sup>6</sup>

## Antimicrobial Susceptibility Testing

Guidance for testing clinical isolates of *S. schleiferi* are found in the Clinical and Laboratory Standards Institute (CLSI) M100 document and by the European Committee on Antimicrobial Susceptibility Testing (EUCAST).<sup>7,9</sup> Breakpoints are available for disk diffusion and dilution methods.<sup>7,9</sup> A CLSI working group evaluated phenotypic methods and breakpoints to determine *mecA*-mediated oxacillin resistance in a set of human and veterinary *S. schleiferi* isolates.<sup>6</sup> The authors showed that the best methods to determine the presence of *mecA*-mediated resistance in *S. schleiferi* was to use oxacillin disk diffusion or minimal inhibitory concentration (MIC) testing applying the *S. pseudintermedius* M100-S28 interpretations.<sup>6,8</sup> Cefoxitin disk diffusion and MIC testing are unreliable for detection of *mecA*-mediated resistance in *S. schleiferi* and should not be used clinically due to high numbers of false-susceptible results.<sup>6,7</sup> EUCAST recommends the use of an oxacillin screen performed by disk diffusion (1µg oxacillin disk) to detect *mecA*-mediated resistance in *S. schleiferi*.<sup>9</sup>

## Key Points

- *Staphylococcus schleiferi* is an increasingly encountered zoonotic pathogen that colonizes the skin and mucosal membranes of small animals, including dogs.
- Due to β-hemolysis on blood agar and variable positivity by coagulase tests (eg, slide, tube and latex agglutination methods), misidentification of *S. schleiferi* as *S. aureus* or other *Staphylococcus* species can occur.
- MALDI-TOF MS is a reliable method to identify and differentiate *S. schleiferi* from other *Staphylococcus* species. Differentiation of staphylococci to the species level is growing increasingly important due to differences in methods and breakpoints applied to detect *mecA*-mediated oxacillin resistance among staphylococci, including *S. schleiferi*.

## References

1. Becker K, Skov RL, Von Eiff C. *Staphylococcus, Micrococcus, and Other Catalase-Positive Cocci*. Ch. 22 in: Carroll KC, Pfaller MA, Landry ML, McAdam AJ, Patel R, Richter SS, Warnock DW (eds), *Manual of Clinical Microbiology*. ASM Press;2019: 367-398.
2. Yarbrough ML, Hamad Y, Burnham CAD, et al. The Brief Case: Bacteremia and Vertebral Osteomyelitis Due to *Staphylococcus schleiferi*. *J Clin Microbiol*. 2017;55:3157–3161.
3. van Griethuysen A, Bes M, Etienne J, et al. International multicenter evaluation of latex agglutination tests for identification of *Staphylococcus aureus*. *J Clin Microbiol*. 2001;.39:86–89. doi:10.1128/JCM.39.1.86-89.2001
4. Delport JA, Peters G, Diagre D, et al. Identification of coagulase-negative staphylococci by the Bruker MALDI-TOF Biotyper compared to the Vitek 2 and MIS gas liquid chromatography. *J Bacteriol Mycol*. 2015;Open Access1(1):00003. doi:10.15406/jbmoa.2015.01.00003.
5. Argemi X, Riegel P, Lavigne T, et al. Implementation of matrix-assisted laser desorption ionization–time of flight mass spectrometry in routine clinical laboratories improves identification of coagulase-negative staphylococci and reveals the pathogenic role of *Staphylococcus lugdunensis*. *J Clin Microbiol*. 2015;53:2030–2036. doi:10.1128/JCM.00177-15.
6. Huse HK, Miller SA, Chandrasekaran S, et al. Evaluation of Oxacillin and Cefoxitin Disk Diffusion and MIC Breakpoints Established by the Clinical and Laboratory Standards Institute for Detection of *mecA*-Mediated Oxacillin Resistance in *Staphylococcus schleiferi*. 2018.
7. Clinical and Laboratory Standards Institute. Performance Standards for Antimicrobial Susceptibility Testing.

Thirtieth Informational Supplement. M100-S30. CLSI; 2020.

8. Clinical and Laboratory Standards Institute. Performance Standards for Antimicrobial Susceptibility Testing. Twenty-Eighth Informational Supplement. M100-S28. CLSI; 2018.
9. European Committee on Antimicrobial Susceptibility Testing, European Society of Clinical Microbiology and Infectious Diseases. [www.eucast.org/clinical\\_breakpoints/](http://www.eucast.org/clinical_breakpoints/)