



COLLEGE of AMERICAN
PATHOLOGISTS

Less is More


Real-World of Digital Pathology and AI

Diana M Lin, MD, FCAP
Dibson D. Gondim, MD, FCAP
Eyas M Hattab MD, MBA, FCAP

November 7, 2023

Diana Murro Lin MD, FCAP

- **Member– Practice Management Committee**
- **PMC Liaison –Digital Content Committee**
- **Associate Professor of Pathology – University of Alabama at Birmingham**
- **Medical Director, AP Division, UAB Highlands Hospital**
- **CAP State Commissioner for Alabama and Division Commissioner for the Florida Panhandle**
- **Board certified AP/CP and Cytopathology**



Eyas M Hattab MD, MBA, FCAP

- **Member, Practice Management Committee, House of Delegates**
- **Professor and Chair of Pathology and Laboratory Medicine at University of Louisville**
- **AP/CP certified with subspecialty fellowships in neuropathology and oncological surgical pathology**
- **Internationally recognized neuropathologist**



Dibson Dibe Gondim, MD, FACP

- **Member of the Artificial Intelligence Committee, CAP**
- **Director of Pathology Informatics and Assistant Professor of Pathology at the University of Louisville**
- **Certified by the American Board of Pathology in:**
 - **Anatomic Pathology**
 - **Neuropathology**
 - **Clinical Informatics**
- **Lead on the large-scale digital pathology and AI Initiative at the University of Louisville**



Disclaimer

The information presented today represents the opinions of the panelists and does not represent the opinion or position of the CAP.

This should not be used as a substitute for professional assistance.

The information in this presentation is provided for educational purposes only and is not legal advice.

Topics for Today's Discussion

- Introduction
- CLIA updates
- Billing
- Real world applications of digital pathology
- Artificial intelligence
- Integrated sign-out dashboard

What is digital pathology?

- Acquiring and interpreting pathology information from a digitized glass slide
 - Diagnosis
 - Consultation
 - Teaching
 - Multi-disciplinary conference presentation
 - And so much more!

Artificial Intelligence (AI)

- Developing algorithms and models using digital analysis and machine learning

Is remote sign-out feasible post COVID?



Yes, you can sign out digital slides without a CLIA certificate for your house or other remote location!



As long as the primary lab meets updated guidelines

Memorandum: "Clinical Laboratory Improvement Amendments of 1988 (CLIA) Post-Public Health Emergency (PHE) Guidance" Center for Clinical Standards and Quality/Quality, Safety & Oversight Group. Centers for Medicare & Medicaid Services. May 11, 2023.

Primary site requirements

- **Has CLIA certification and complies with all relevant Federal laws.**
- **Is certified to perform all the of work performed at the remote site and the laboratory director is responsible for all testing.**
- **Retains all documentation, including tests performed remotely and list of staff working remotely.**
- **Indicates remote site location on the reports (coding system ok).**

Memorandum: "Clinical Laboratory Improvement Amendments of 1988 (CLIA) Post-Public Health Emergency (PHE) Guidance" Center for Clinical Standards and Quality/Quality, Safety & Oversight Group. Centers for Medicare & Medicaid Services. May 11, 2023.

Yes, you can code!

- **Category III digital CPT codes now available**
- **Must use whole-slide imaging**
- **No codes for**
 - Education/research
 - Tumor boards/clinical conferences
 - Archiving
 - Validation

Diagnosis, Special Stains and IHC

Digital Code	In conjunction with
+0751T	88302
+0752T	88304
+0753T	88305
+0754T	88307
+0755T	88309
+0756T	88312
+0757T	88313
+0758T	88314
+0759T	88319
+0760T	88342
+0761T	88341
+0762T	88344
+0763T	88360

Coming soon on Jan. 1- cytology, expanded surgical and hematology codes

Digital Code	In conjunction with	Digital Code	In conjunction with	Digital Code	In conjunction with
+0827T	88104	+0838T	88321	+0849T	88364
+0828T	88106	+0839T	88323	+0850T	88366
+0829T	88108	+0840T	88325	+0851T	88368
+0830T	88112	+0841T	88331	+0852T	88369
+0831T	88141	+0842T	88332	+0853T	88377
+0832T	88160	+0843T	88333	+0854T	85060
+0833T	88161	+0844T	88334	+0855T	85097
+0834T	88162	+0845T	88346	+0856T	88348
+0835T	88172	+0846T	88350		
+0836T	88177	+0847T	88363		
+0837T	88173	+0848T	88365		

Digital Pathology: Real-world applications


- Not one size fits all
- Many different applications
 - Specific operational / logistical function
 - Specific mission (education, research, etc)
 - Improved staffing / access
 - Grow business (referrals, consultations)
- Means to an end! (AI deployment)

Inhouse consultation

SLIDE TRAY Panel: Slide 1


A: BRAIN 1

HE




I have an interesting case and would appreciate your expert opinion.

0.4x



Of course!




Thank you. Here is the link:
<https://digitalslide/case>



I see. Is there something specific that concerns you?



Actually, yes. I can set the concerning area as the start view. Here:
<https://digitalslide/case/zoomLevel>



Great, I will take a look...

Inhouse consultation



I have an interesting case and would appreciate your expert opinion.

Of course!



Thank you. Here is the link:
<https://digitalslide/case>

I see. Is there something specific that concerns you?



Actually, yes. I can set the concerning area as the start view. Here:
<https://digitalslide/case/zoomLevel>

Great, I will take a look...



- Easily communicate in real time!
- No shuffling slides between locations!
- Better turnaround time!

External digital consultation / point of care

- Using digital pathology for external consultation
 - ~~Cumbersome physical workflow~~
 - ~~Extra paperwork~~
 - ~~Slides packaged / sent~~
 - ~~Shipping costs and labor add up over time~~
 - ~~Suboptimal turnaround time~~
 - ~~Original slides lost, broken, or never returned~~

Conferences and Tumor Board

- No more pulling slides!
- No more photographing the slides!
- No more PowerPoint presentations!



Consensus Conference



Image from zeiss.com



Consensus Conference




- Better attendance!
- No more crowding!
- Sign in from anywhere!



“Tagging” cases

- **Advanced cataloging / organizing**
- **Automated workflow processes**
- **Easy data sharing and collaboration**
- **A tag can be any label! Customizable**

“Tagging” cases



Tagging through the digital slide viewer

“Tagging” cases

Pathology Report Internal Comments

ssss/+40+SP23000XXX001001001/pituitary fossa/pitnet, corticotroph type, sparsely granulated, silent

accession	pathologist	slide_site	slide_comment	slide_url
03-SP-23-XXX	HATTAB, EYAS, MD-PAT	cornea	fungal keratitis	https://DigitalSlideViewer/launchCase
03-SP-23-XXX2	HATTAB, EYAS, MD-PAT	cornea	chronic corneal edema with bullous keratopathy	https://DigitalSlideViewer/launchCase
40-SP-23-XX	HATTAB, EYAS, MD-PAT	clivus	rathke cleft cyst	https://DigitalSlideViewer/launchCase
40-SP-23-XX	HATTAB, EYAS, MD-PAT	sphenoid mass	chordoma	https://DigitalSlideViewer/launchCase
03-SP-23-X3	HATTAB, EYAS, MD-PAT	eyelid	merkel cell carcinoma	https://DigitalSlideViewer/launchCase
40-SP-23-XXX	HATTAB, EYAS, MD-PAT	pituitary fossa	pitnet, corticotroph type, sparsely granulated, silent	https://DigitalSlideViewer/launchCase
40-SP-22-X1	HATTAB, EYAS, MD-PAT	spine	myxopapillary ependymoma	https://DigitalSlideViewer/launchCase

Tagging through the pathology report to send to a common database

“Tagging” cases

to photograph

potential GSU

tumor board

resident lecture

case of interest

UPDATED INTEGRATED DIAGNOSIS

FINAL INTEGRATED DIAGNOSIS:

Solitary fibrous tumor, CNS WHO grade 2

HISTOLOGIC DIAGNOSIS:

Spindle cell tumor

RELEVANT ANCILLARY INFORMATION:




STAT6 (IHC): Positive

Consensus Methylation Profiling Class: Solitary fibrous tumor, CNS WHO grade 2 (see Comment)

CNS WHO GRADE:

WHO grade 2

Case View

		
<p>Block: A - 001 - 001</p> <p>Block-Description: Left Lateral Orbital Mass</p> <p>Case: Surgical Pathology</p>	<p>Block: A - 001 - 002</p> <p>Block-Description: Left Lateral Orbital Mass</p> <p>Case: Surgical Pathology</p>	<p>Block: A - 002 - 001</p> <p>Block-Description: Left Lateral Orbital Mass</p> <p>Case: Surgical Pathology</p>


Tagging through other sources

DP Implementation is not without challenges

- **Cost/Business case (initial investment and ongoing maintenance)**
- **Regulatory issues**
- **Buy in (C-suite, pathologists, clinicians, staff, etc)**
- **Requires significant workflow adjustments**
- **Requires IT support, clinical information management, etc**
- **Vulnerable to interruptions**
- **Scanning/viewing limitations**
 - Cytology, bone marrows, H. pylori, mitoses, polarization, etc
- **Requires a QA program**
 - Failure rate, re-scanning

Artificial Intelligence

- Basic requirements
- Common AI-based systems for histopathology
- Ensuring trustworthiness in AI pathology
- Maximizing productivity
- AI + DP interoperability



Pathologists: The pillars of precision medicine

- **Masters of doing more with less**
- **Efficiency masters**
 - Unparalleled skill in maximizing precision with fixed resources
- **Precision medicine champions**
 - Crucial to the seamless integration of advanced diagnostics into patient care
- **Escalating duties**
 - Steering through intricate molecular diagnostics landscape
 - Absorbing rapid growth in medical knowledge
 - Adapting to a tight web of regulatory requirements
- **Sustainability check**
 - Evaluating the potential for overload
 - Necessity of innovative tools to support pathologists' expanding roles



The role of AI in pathology

- **AI as a tool**
 - Requires careful exploration and mastery
 - Still in early stages—potential unknown risks
 - Our duty to understand and mitigate risks
 - Robust DP Infrastructure is critical
- **Implementation challenges**
 - Substantial effort for AI system deployment
 - Not universally established infrastructure
- **Productivity and AI**
 - AI deployment alone ≠ guaranteed productivity gains
 - True potential lies in AI + knowledge + optimal implementation

Common AI-based systems for histopathology



- #1** Classification & detection
- #2** Digital biomarker and virtual staining
- #3** Image-based predictive/prognostic model

(Note: This list is not exhaustive)

#1 Classification & detection systems

- **Tumor detection, quantification, and grading**
 - Advanced algorithms for tumor analysis
 - Example: solutions for prostate and breast biopsies
- **Metastasis detection**
 - Automated detection of lymph node metastases
- **Mitotic activity assessment**
 - AI-driven counting of mitotic figures
- **Result verification**
 - Tools for pathologists to confirm AI findings
 - Visualization techniques
 - Use of masks and heatmaps for AI-histology correlation

#1 Example: AI for cancer identification, quantification, grading, and PNI detection



20x

AI


AI PAIGE PROSTATE

RUO

Focus of Interest

Invasive Carcinoma	
Invasive Carcinoma	62.71 %
Length	8.65 mm
Pattern 4	100.00 %
Gleason Score	4 + 4 = 8
Perineural Invasion	

#1 Example: AI for kidney tumor classification with heatmaps




More examples can be found here: <http://aistain.com>

Gondim DD, Al-Obaidy KI, Idrees MT, Eble JN, Cheng L. Artificial intelligence-based multi-class histopathologic classification of kidney neoplasms. *J Pathol Inform.* 2023 Feb 16;14:100299. doi: 10.1016/j.jpi.2023.100299. PMID: 36915914; PMCID: PMC10006494.

#2 Digital biomarker and virtual staining

- **Enhanced turnaround time**
 - Accelerate result availability, bypassing the wait for traditional staining or other ancillary studies.
- **Resource efficiency**
 - Decreases dependence on physical instruments, supplies and specialized personnel
- **Consideration of risks**
 - Validation could be performed based on gold standard
 - Challenge prospective monitoring
 - Limited literature about real world deployment

#2 Example: Virtual trichrome for liver



Naglah A, Khalifa F, El-Baz A, Gondim D. Conditional GANs based system for fibrosis detection and quantification in hematoxylin and eosin whole slide images. Medical Image Analysis. 2022 Oct 1;81:102537.

*Systems and methods for digital transformation of medical images and fibrosis detection
Provisional patent, US Patent App. 17/845,880
El-Baz, A.S., Gondim, D., Naglah, A., and Khalifa, F.*

#3 Image-based predictive/prognostic models

- **Optimal treatment guidance**
 - Utilizing AI to recommend the most effective treatment options tailored to individual patient profiles
- **Prognosis group stratification**
 - AI algorithms assist in categorizing patients into prognosis groups, enabling personalized care plans
- **Case study: ArteraAI prostate test**
 - An example of AI's application in providing prognostic insights for prostate cancer management
- **Verification challenges**
 - The current limitations in secondary methods for independent result confirmation


Ensuring trustworthiness in AI pathology

- **Critical role of validation**
 - Validation protocols are paramount to establish the reliability of AI results
- **Pathologist oversight (human-in-the-loop)**
 - Pathologist in the decision-making loop to adjudicate AI findings
 - Pathologist oversight not possible for biomarkers, virtual staining, and predictive/prognostic assessments
- **Adjudication workflow integration**
 - Incorporating a structured workflow for pathologists to review and adjust AI outcomes
- **Risk mitigation**
 - Addressing the potential for false negatives or positives, which could inadvertently increase the pathologist's workload

DP + AI Interoperability

- **Full-spectrum integration (ideal)**
 - EHR + LIS + DP + AI
 - Lab-centric control (LIS + DP + AI)
- **Efficiency in case management**
 - Eliminates the need for manual retrieval across different applications
 - Allows order placement from image management systems (IMS) to laboratory information system (LIS)
 - Easy export of images from digital pathology (DP) to LIS
- **Seamless asset tracking**
 - Simplifies the monitoring of asset status, improving turnaround times and resource allocation


DP + AI Interoperability



Standard-based interfaces


Potential to plug and play

(NOT A REALITY IN THE MARKET)




Custom interfaces

(Not scalable – Not sustainable)




Custom interfaces

(Not scalable – Not sustainable)




Custom interfaces

(Not scalable – Not sustainable)



Example: LIS + DP + AI Interoperability



Conclusion & take-home messages

- DP and AI offer the potential for enhanced accuracy and productivity
- High-quality histology and lab operational efficiency are paramount
- Leveraging AI's full potential
 - Expertise, exploration, validation, and interoperability
- Adjudication risk
 - **Suboptimal accuracy = workflow complexities and higher workloads**
- Interoperability risk
 - **Inadequate interoperability = workflow complexities and higher workloads**



COLLEGE of AMERICAN
PATHOLOGISTS

Membership


Did you find this information useful?

This program was funded by your CAP membership. Please be sure to keep your membership current so we can continue to bring timely and relevant resources like this to you.

Visit cap.org to renew your membership or email membership@cap.org.

Pathology Business Fundamentals

Essential online courses to help grow your management skills to lead your practice



1. Relative Value Units (RVU's)—Understanding the Basics
2. How Pathology Practices Get Paid
3. Revenue Cycle Management
4. Analysis and Interpretation of Billing Reports
5. Basic Practice Cost Analysis
6. Capacity Management and Workflow Analysis
7. Basic Contracting and Fee Analysis
8. Basic Budget Development

Learn more
and register



Additional Resources

Practice Management

- <https://www.cap.org/member-resources/practice-management>

Practice Management Articles

- <https://www.cap.org/member-resources/articles/category/practice-management>

We value your feedback!



If after attending this discussion and later you applied any of what you learned to your practice, please share your feedback of how it worked for your practice at <https://www.cap.org/member-resources/practice-management/practice-management-inquiry-form> .



Watch for the session evaluation form. Your feedback is important!