


1

Today's Presenter

Earle Collum, MD, FCAP

- Earle Collum MD, FCAP, is the Division Medical Director for Mid-America as well as the site medical Director for the Phoenix and Denver regional labs at LabCorp. Previously, he was the Medical Director of St Joseph's Hospital and Medical Center's Department of Pathology and Laboratory Medicine in Phoenix. He was also President of the Medical Staff. Dr. Collum is Board Certified in Anatomic and Clinical Pathology. He is an active member of the CAP and has performed over 200 inspections. Currently he is the chair of the CAP's Complaints and Investigations Committee. He was recognized with a CAP Lifetime Achievement award in 2017.




© College of American Pathologists 8 June 2020 2

2

Today's Presenter

Richard M. Scanlan, MD, FCAP

- Richard M. Scanlan, MD, FCAP, is vice chair of clinical pathology, professor of pathology, and director of transfusion medicine at Oregon Health Sciences University (OHSU) in Portland, Oregon. He is board certified in anatomic and clinical pathology. He is actively involved in the Portland pathology community and serves as laboratory director for Planned Parenthood of Oregon. Dr. Scanlan is a CAP governor and has been active in the CAP's Laboratory Accreditation Program since 1986. Currently he is the chair of the CAP Council on Laboratory Accreditation.




© College of American Pathologists 8 June 2020 3

3

Objectives

- Common COVID-19 Terminology
- Current COVID-19 Tests
- Verification/validation of COVID-19 testing
- Quality monitoring and personnel qualifications
- Biosafety for COVID-19
- General COVID-19 testing questions



4

Terminology

5

Emergency Use Authorization (EUA)

- What is an EUA?
 - Legal mechanism used by the FDA
 - Allows use of unapproved medical products during a public health emergency
 - Used when there are no adequate, approved, and/or available alternatives

6

Emergency Use Authorization (EUA) – cont'd.

- How is an EUA different than FDA-clearance or approval?
- What about laboratories outside of the US?

© College of American Pathologists

7

7

Laboratory-Developed Test (LDT)

- What is an LDT?
 - Test is performed by the clinical laboratory in which it was developed wholly or in part AND
 - Test is neither FDA cleared/approved



© College of American Pathologists

8

8

Laboratory-Developed Test (LDT) – cont'd

- Is an EUA considered an LDT?
- Can my laboratory develop an LDT for COVID-19 testing?


© College of American Pathologists

9

9

Test Complexity

- What is the complexity of tests used for COVID-19?
 - Prior to obtaining EUA
 - Without EUA (following FDA policy)
 - With EUA



© College of American Pathologists 10

10

Test Complexity cont.

- Are there any waived tests?

* A complete listing of EUA tests and the complexity of each is found on the [FDA website](#).

© College of American Pathologists 11

11

Types of Testing

© College of American Pathologists 12

12

Types of COVID-19 Testing

- What types of COVID-19 tests are available?
 - Molecular diagnostic tests
 - Serological tests
 - Antigen detection tests

* Helpful references

[The College of American Pathologists: Clarifying Tests for COVID-19](#)


[FDA Policy for Coronavirus Disease-19 Tests During the Public Health Emergency](#)

© College of American Pathologists

13

13

Serology Testing




© College of American Pathologists

14

14


Molecular Testing



© College of American Pathologists

15


15



16



17



18

Modifications to EUAs

What does my laboratory do if there is a shortage of essential component needed to perform the testing?

- Supply shortages have required alternatives
 - Swabs
 - Transport media
 - Testing materials
 - Equipment


© College of American Pathologists

19

19

Modifications to EUAs cont.

- FDA provides recommendations on alternatives in their [FAQs](#)
- Laboratories must validate other alternatives (bridging study)



© College of American Pathologists

20

Quality Monitoring

© College of American Pathologists

21

21

Proficiency Testing (PT) – COM.01500

- Does the CAP have PT available for COVID-19 testing?
 - Molecular-based testing
 - Serologic testing
 - Antigen testing

© College of American Pathologists

22

22

Proficiency Testing (PT) – COM.01500

- Perform alternative performance assessment if the laboratory is not performing PT

- * Order CAP PT products in the online store or by contacting the CAP at 800-323-4040

Home > Laboratory Development > Proficiency Testing > SARS-CoV-2 Proficiency Testing Programs

SARS-CoV-2 Proficiency Testing Programs

As laboratories increase their testing capabilities for COVID-19, we have implemented the following PT programs to support your efforts.



Enroll Today

Click your order today to receive material immediately.

Order SARS-CoV-2 Serologic (COVID)

Order SARS-CoV-2 Molecular (COVID)

© College of American Pathologists

23

23

Quality Control and IQCP

- What type of QC does my laboratory need to do for COVID-19 testing?
- Does my laboratory need to implement an IQCP for COVID-19 testing?

© College of American Pathologists

24

24

Quality Control and IQCP cont.

- Waived tests (authorized for the patient care setting)
 - Follow manufacturer's instructions for internal/external QC – IQCP does not apply

© College of American Pathologists

25

25

Quality Control and IQCP cont.

- Moderate and high complexity tests
 - If manufacturer instructions define reduced external QC frequency, follow instructions – no IQCP required
 - If manufacturer's instructions do not define reduced external QC frequency
 - Perform positive and negative external QC each day of use **OR**
 - Develop and implement an IQCP



© College of American Pathologists

26

26

Nucleic Acid Contamination and Molecular Testing

- Are there other types of quality monitoring my laboratory should do for COVID-19 testing?



© College of American Pathologists

27

Nucleic Acid Contamination and Molecular Testing cont.

- Laboratories must monitor and minimize contamination to prevent false positive results
 - Perform wipe (environmental) tests
 - Monitor statistics
 - Investigate any physician inquiries
 - Use process controls to minimize risk of contamination
 - during collection, transport, processing and storage
 - Following instructions for handling waste (used test cartridges)


© College of American Pathologists

28

28

Instrument/Equipment Maintenance

- What do we do if scheduled equipment maintenance was not done as scheduled due to the COVID-19 pandemic and vendors not being able to come onsite?



© College of American Pathologists

29

29


Instrument/Equipment Maintenance – cont'd

- The CAP recommends the following actions:
 - Create a record to document and track the issue
 - Identify reasonable alternatives to complete the maintenance or perform an assessment of the affected equipment to determine continued use
 - Monitor continued performance through QC/function checks
 - Work with vendors to complete necessary maintenance
 - Retain records of the issue and resolution for at least two years


© College of American Pathologists

30


30




31




32




33



34



35



36

Personal Protective Equipment (PPEs) cont.

- What precautions should be used when handling/testing specimens from suspected or confirmed COVID-19 patients?



© College of American Pathologists

37

37

Personal Protective Equipment (PPEs) cont.

- For handling/testing, CDC recommends
 - Standard precautions including gloves, lab coat/gown, eye protection
 - For procedures with a **high likelihood to generate droplets or aerosols**
 - Surgical mask/face shield
 - Other physical barriers
 - Centrifuge covers/caps
 - Biologic safety cabinets

* Visit the [CDC website](https://www.cdc.gov) for the most current guidance

© College of American Pathologists

38

38

Q & A'S on COVID-19 Testing

© College of American Pathologists

39

39

Recommended Testing

- Does the CAP recommend any specific COVID-19 testing for my laboratory setting?
 - 1) test volumes
 - 2) turn around times needed
 - 3) current testing platforms in use within your laboratory
 - 4) testing personnel qualifications
 - 5) PPE and biosafety measures available for staff
 - 6) availability of supplies to begin and maintain patient testing

© College of American Pathologists

40

40

Pooling of Samples

- Can my laboratory perform testing on pools of specimens for COVID-19 molecular-based testing?

* Contact the FDA Division of Microbiology for additional information on regulatory approval for modifying an EUA

- (301) 348-1778
- CDRH-EUA-Templates@fda.hhs.gov

© College of American Pathologists

41

41

Reporting of Results for Multiple Platforms

- My laboratory is using three different instrument platforms for COVID-19 molecular testing and each uses different reporting language. We use the same reportable result for all three platforms?



© College of American Pathologists

42

Specimen Collection Devices

- We have two different instruments for COVID-19 EUA molecular testing. The manufacturer's instructions for one includes the use of nasopharyngeal swabs and nasal swabs. The other only includes nasopharyngeal swabs. Can we validate the second test to allow the use of nasal swabs on both instruments?

© College of American Pathologists

43

43

Transporting and Shipping Specimens

- How should suspected or confirmed COVID-19 specimens be transported within a facility?
- What are the requirements for packaging and shipping specimens from suspected or confirmed COVID-19 patients?

© College of American Pathologists

44

44

Temporary Suspension of Testing

- Due to the COVID-19 pandemic, my laboratory temporarily suspended some of its testing. What do we need to do restart it?

© College of American Pathologists

45

45

Disclaimers in reporting

- Are specific disclaimers required on reports for EUA COVID-19 tests?

© College of American Pathologists

46

46

Competency Assessment

- With the current challenges due to the public health crisis, is any leniency being given to performing competency assessment?
- How can we perform direct observation with social distancing?



© College of American Pathologists

47

47

New Reagent Lot Testing

- My laboratory has two identical instruments for COVID-19 testing. Do we need to do new reagent lot testing on both?
- If my laboratory has the same instrument as our satellite laboratory, can we QC the new lot and send it to the other site to use also?

© College of American Pathologists

48

48

Problems with EUA Tests

- What actions should my laboratory take if we notice problems with the EUA test it is using, such as complaints from physicians about false positives or negatives?

© College of American Pathologists

49

49

Laboratory Director Sign Off

- My laboratory director has been unable to come into the laboratory to sign off on proficiency testing, test validation/verification studies, and other items. Can we obtain the signature electronically?



© College of American Pathologists

50

50

Missed/Delayed Duties

- Because some staff were furloughed or unable to come to work, there are some instances where my laboratory missed performing required duties (eg, instrument maintenance, PT, method comparison) or performed them later than required. Will my laboratory be cited during our next inspection?

© College of American Pathologists

51

51

Inspection Delays

- My laboratory's CAP inspection window is due to open in August 2020, should we expect any delays?
- My laboratory's anniversary date was in May 2020 and we have not had our inspection yet. Are we still considered CAP-accredited?



© College of American Pathologists

52

Direct to Consumer Testing (DTC)

Can we perform DTC for COVID-19 testing?

- State by State

© College of American Pathologists

53

53

Workplace or Corporate Testing


Can we test for a corporation and report results back to its HR department?



© College of American Pathologists

54

54



55



56



57
